

**Model 276  
Photomultiplier Base  
Operating and Service Manual**

# **Advanced Measurement Technology, Inc.**

a/k/a/ ORTEC<sup>®</sup>, a subsidiary of AMETEK<sup>®</sup>, Inc.

## **WARRANTY**

ORTEC\* warrants that the items will be delivered free from defects in material or workmanship. ORTEC makes no other warranties, express or implied, and specifically NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

ORTEC's exclusive liability is limited to repairing or replacing at ORTEC's option, items found by ORTEC to be defective in workmanship or materials within one year from the date of delivery. ORTEC's liability on any claim of any kind, including negligence, loss, or damages arising out of, connected with, or from the performance or breach thereof, or from the manufacture, sale, delivery, resale, repair, or use of any item or services covered by this agreement or purchase order, shall in no case exceed the price allocable to the item or service furnished or any part thereof that gives rise to the claim. In the event ORTEC fails to manufacture or deliver items called for in this agreement or purchase order, ORTEC's exclusive liability and buyer's exclusive remedy shall be release of the buyer from the obligation to pay the purchase price. In no event shall ORTEC be liable for special or consequential damages.

### **Quality Control**

Before being approved for shipment, each ORTEC instrument must pass a stringent set of quality control tests designed to expose any flaws in materials or workmanship. Permanent records of these tests are maintained for use in warranty repair and as a source of statistical information for design improvements.

### **Repair Service**

If it becomes necessary to return this instrument for repair, it is essential that Customer Services be contacted in advance of its return so that a Return Authorization Number can be assigned to the unit. Also, ORTEC must be informed, either in writing, by telephone [(865) 482-4411] or by facsimile transmission [(865) 483-2133], of the nature of the fault of the instrument being returned and of the model, serial, and revision ("Rev" on rear panel) numbers. Failure to do so may cause unnecessary delays in getting the unit repaired. The ORTEC standard procedure requires that instruments returned for repair pass the same quality control tests that are used for new-production instruments. Instruments that are returned should be packed so that they will withstand normal transit handling and must be shipped PREPAID via Air Parcel Post or United Parcel Service to the designated ORTEC repair center. The address label and the package should include the Return Authorization Number assigned. Instruments being returned that are damaged in transit due to inadequate packing will be repaired at the sender's expense, and it will be the sender's responsibility to make claim with the shipper. Instruments not in warranty should follow the same procedure and ORTEC will provide a quotation.

### **Damage in Transit**

Shipments should be examined immediately upon receipt for evidence of external or concealed damage. The carrier making delivery should be notified immediately of any such damage, since the carrier is normally liable for damage in shipment. Packing materials, waybills, and other such documentation should be preserved in order to establish claims. After such notification to the carrier, please notify ORTEC of the circumstances so that assistance can be provided in making damage claims and in providing replacement equipment, if necessary.

---

Copyright © 2005, Advanced Measurement Technology, Inc. All rights reserved.

\*ORTEC<sup>®</sup> is a registered trademark of Advanced Measurement Technology, Inc. All other trademarks used herein are the property of their respective owners.

## CONTENTS

WARRANTY .....	ii
SAFETY INSTRUCTIONS AND SYMBOLS .....	iv
SAFETY WARNINGS AND CLEANING INSTRUCTIONS .....	v
1. DESCRIPTION .....	1
2. SPECIFICATIONS .....	1
2.1. PHOTOMULTIPLIER TUBE BASE .....	1
2.1.1. PERFORMANCE .....	1
2.1.2. CONTROL .....	1
2.1.3. INPUT .....	1
2.1.4. OUTPUTS .....	1
2.2. PREAMPLIFIER .....	1
2.2.1. PERFORMANCE .....	1
2.2.2. OUTPUT .....	2
2.2.3. INPUTS .....	2
2.3. ELECTRICAL AND MECHANICAL .....	2
3. INSTALLATION .....	2
3.1. DETECTOR MOUNTING .....	2
3.2. SYSTEM CONNECTION .....	2
3.3. INITIAL ADJUSTMENTS .....	2
4. OPERATION .....	3
4.1. CALCULATION OF RESPONSE OF SCINTILLATOR/PHOTOMULTIPLIER .....	3
4.2. TIMING APPLICATIONS .....	5
4.2.1. TYPICAL FAST-SLOW COINCIDENCE SYSTEM .....	5
4.2.2. GAMMA-GAMMA COINCIDENCE SYSTEM FOR THE HIGH-PURITY GERMANIUM (HPGe) DETECTOR .....	6
4.3. SCINTILLATION SPECTROSCOPY .....	6
5. MAINTENANCE .....	7
5.1. FACTORY REPAIR .....	8
6. BIBLIOGRAPHY .....	9

## SAFETY INSTRUCTIONS AND SYMBOLS

This manual contains up to three levels of safety instructions that must be observed in order to avoid personal injury and/or damage to equipment or other property. These are:

**DANGER** Indicates a hazard that could result in death or serious bodily harm if the safety instruction is not observed.

**WARNING** Indicates a hazard that could result in bodily harm if the safety instruction is not observed.

**CAUTION** Indicates a hazard that could result in property damage if the safety instruction is not observed.

Please read all safety instructions carefully and make sure you understand them fully before attempting to use this product.

In addition, the following symbol may appear on the product:



**ATTENTION – Refer to Manual**



**DANGER – High Voltage**

Please read all safety instructions carefully and make sure you understand them fully before attempting to use this product.

## SAFETY WARNINGS AND CLEANING INSTRUCTIONS

**DANGER** Opening the cover of this instrument is likely to expose dangerous voltages. Disconnect the instrument from all voltage sources while it is being opened.

**WARNING** Using this instrument in a manner not specified by the manufacturer may impair the protection provided by the instrument.



### Cleaning Instructions

To clean the instrument exterior:

- Unplug the instrument from the ac power supply. Remove loose dust on the outside of the instrument with a lint-free cloth.
- Remove remaining dirt with a lint-free cloth dampened in a general-purpose detergent and water solution. Do not use abrasive cleaners.

**CAUTION** To prevent moisture inside of the instrument during external cleaning, use only enough liquid to dampen the cloth or applicator.

- Allow the instrument to dry completely before reconnecting it to the power source.



# ORTEC MODEL 276 PHOTOMULTIPLIER BASE

## 1. DESCRIPTION

The ORTEC 276 Photomultiplier Base is designed to provide a near-optimum voltage distribution for essentially all 10-stage photomultiplier tubes that fit the standard 14-pin tube socket. The unit has a self contained low-noise integrating preamplifier that provides gain and line drive capability for the output signal from the last dynode. The preamplifier is dc-coupled to simplify pole-zero cancellation in the main amplifier. The anode signal is capacitively coupled to an output BNC connector for use in coincidence timing or linear pulse-height analysis. A TEST input is provided to allow pulser calibration and testing. A FOCUS control on the rear panel is for optimizing the photocathode first-dynode electrical field. Some of the photomultiplier tubes with which this PM base is compatible are:

- **RCA:** 4518, 5819, 6217, 6342A, 6655A, 7326, 8053, 8054, and 8055
- **Philips:** XP-1000 to 1005 and XP-1031 through 1033
- **EMI:** 9536, 9578, 9579, and 9708 Series
- **CBS:** 7817, 7818, 7819, and CL-1004 to 1012
- **Dumont:** 6292, 6363, and 6364

The unit is also compatible with other 10-stage tubes not listed above. Compatibility can be determined by comparison with those listed.

## 2. SPECIFICATIONS

### 2.1. PHOTOMULTIPLIER TUBE BASE

#### 2.1.1. PERFORMANCE

**Design** Resistor divider bleeder connected to 14-pin standard PM tube base for 10-stage tubes. Total resistance 1.49 M $\Omega$ , with bleeder current of 1.33 mA when maximum high voltage of 2000 V is applied.

**Voltage Distribution Coefficient** Linear to all stages with focus adjustment on the grid.

#### 2.1.2. CONTROL

**FOCUS** Single-turn locking potentiometer on panel for external adjustment of PM tube grid potential.

#### 2.1.3. INPUT

**HIGH VOLTAGE** Type SHV connector on panel accepts operating voltage for distribution to PM tube; +2000 V maximum.

#### 2.1.4. OUTPUTS

**ANODE** Type BNC connector on panel furnishes PM tube anode negative output pulse for use primarily as a timing signal;  $Z_o = 1$  k $\Omega$  ac-coupled.

**LAST ANODE** Connected internally to preamplifier input.

### 2.2. PREAMPLIFIER

#### 2.2.1. PERFORMANCE

**Input** Directly from last dynode of photomultiplier tube.

**Conversion Gain** Nominally 5  $\mu$ V/eV using a 2  $\times$  2-in. NaI(Tl) crystal and a photomultiplier tube gain of  $10^6$ .

**Rise Time** <100 ns for a fast rise time test pulse.

**Fall Time Constant** 50  $\mu$ s.

**Integral Nonlinearity** < $\pm$ 0.02%, 0 to +10 V.

**Output Polarity** Positive.

**Output Noise** <50  $\mu$ V rms.

**Saturation Level** +10 V into open circuit; +5 V into 93  $\Omega$  load.

**Gain Temperature Coefficient** < $\pm$ 0.005%/°C, 0 to 50°C.

### 2.2.2. OUTPUT

**PREAMP** Type BNC connector on panel furnishes preamplifier positive output pulse for use primarily as an energy analysis signal.  $Z_o = 93 \Omega$  dc-coupled.

### 2.2.3. INPUTS

**Power Cable** 3-m (10-ft) captive cable terminates in Amphenol 17-20090 connector to accept operating power for the preamplifier circuits. Compatible with all ORTEC main amplifiers and with an ORTEC 114 Preamplifier Power Supply.

**TEST** Type BNC connector on panel accepts test pulses from an ORTEC Pulse Generator for testing and calibration. Nominal 100- $\Omega$  charge termination built into preamplifier circuit.

## 2.3. ELECTRICAL AND MECHANICAL

### Power Requirements

**Preamplifier** +24 V, 16 mA; -24 V, 16 mA.

**PM Tube Base** +2000 V maximum (use rate voltage for the PM tube that is installed).

### Weight

**Net** 0.65 kg (1.5 lb).

**Shipping** 1.3 kg (3.0 lb).

**Dimensions** 5.58 cm (2.2 in.) diam  $\times$  10.1 cm (4.0 in.) long plus 3-m (10-ft) captive power cable.

## 3. INSTALLATION

### 3.1. DETECTOR MOUNTING

Normally with 10-stage photomultipliers the amount of signal current through the tube is small; therefore, it is quite practical to run positive high voltage on the tube with the cathode at ground potential. This means that when the detector is mounted to the photocathode, very little attention need be paid to the probability of leakage currents being created across the glass envelope, since the voltage potential should be zero. In many cases where this PM base is used, the detector is mounted to the photomultiplier in an integral package.

When a scintillator is mounted on a photomultiplier, the combination must be made light-tight to ensure low noise and to eliminate the possibility of photomultiplier damage. The most common practice to make the detector/photomultiplier light-tight is to carefully wrap it with a good-quality black electrical tape. Another is to use an aluminum can cover and to tape it at the entrance.

### 3.2. SYSTEM CONNECTION

If it is desired to test and calibrate the system, connect a pulser to the test input. A 1-V test pulse should yield a 1-V pulse on the preamplifier output. The test input is internally terminated in 93  $\Omega$ . Connect a high voltage power supply to the high voltage input connector and **ensure positive polarity** of the high voltage supply. Connect the preamplifier signal to the energy analysis system

and connect the anode signal to either a timing or energy analysis system. Use high-quality coaxial cable and terminate the anode signal cable in its characteristic cable impedance.

### 3.3. INITIAL ADJUSTMENTS

The only adjustment on the PM Base is the FOCUS control, which optimizes the photocathode/first-dynode electrical field.

1. Plug a photomultiplier detector into the unit and ensure it is light-tight. Connect the positive high voltage supply.
2. Place a radiation source in the vicinity of the scintillation detector.
3. Observe the preamplifier output on a sensitive range of an oscilloscope and slowly increase the high voltage.
4. If there is a very large amount of unipolar grass (circuit noise), there is possibly a light leak. Cover the PM tube with a black cloth to see whether the grass diminishes; if it does, the photomultiplier was seeing the ambient light and should be rechecked for proper covering to eliminate the light leak.
5. With the high voltage at the desired level, observe the pulses which are radiation-induced scintillations and adjust the FOCUS control to obtain maximum pulse amplitude. This is the proper setting; lock the adjustment with the potentiometer lock nut.

## 4. OPERATION

Once the steps outlined in Section 3 of this manual are performed, the unit is ready for use. High voltage may be applied and adjusted for the appropriate gain associated with the specific experiment. The gain will vary by a factor of approximately 2 for each high-voltage change of 100 V. NOTE: It is advisable to operate the high voltage at the minimum practical value when the high count rates are to be experienced, since count rate tolerance is a direct function of the photomultiplier gain.

### 4.1. CALCULATION OF RESPONSE OF SCINTILLATOR/PHOTOMULTIPLIER

Table 4.1 lists the decay constants of some of the more common scintillators. The first three scintillators are crystals. Naton-136, Pilot B, and NE-102 are plastics, and the last two are liquid scintillators. The decay time  $\tau_1$  is responsible for a finite rise time on the leading edge of L(t) (refs. 12, 21, 22);  $\tau_2$  is the fast decay component which is not noticeable at the output of the photomultiplier;  $\tau_3$  and  $\tau_4$  are the slow components (important for n- $\gamma$  discrimination with NE-213, NE-218, and Stilbene). Where measured values were not given, the letters **N.G.** have been entered. The parameter P is the number of photo electrons released at the photocathode per unit energy. This figure is affected by the efficiency and spectral response of the photocathode (refs. 22, 23, 26, 27) and hence is somewhat characteristic of the photomultiplier used. However, it provides a reasonably good guide for comparing the light output of scintillators. In Table 4.1, P is listed as a fraction of the value for anthracene;  $P_{(\text{anthracene})}$  is  $\sim 700$  eV/photoelectron for S-11 photocathode material.

Table 4.1. Scintillator Decay Times.

Scintillator	Decay Time Constants (ns)				P/P <sub>(anthracene)</sub>
	$\tau_1$	$\tau_2$	$\tau_3$	$\tau_4$	
Nal(Tl)		250	(24)*		2.40
Anthracene (25)*	N.G.	29.3			1.00
Stilbene (23)*	0.1	4.05	33	270	0.64
Naton-136 (22)*	0.4	1.6	10		0.37
Pilot B (22)*	0.4	1.6	9.5		0.41
NE-102 (25)*	N.G.	2.51	N.G.		0.65
NE-213 (23)*	1.66	3.16	32.3	270	0.47
NE-218 (23)*	1.76	3.58	36.5	288	0.56

\*See Bibliography

The thickness of the scintillator is frequently chosen according to the required stopping power. The flight time of the radiation (gamma rays or neutrons) across the scintillator normally becomes a limitation on the time resolution as the thickness is increased. Usually detection efficiency must be compromised for good time resolution.

The scintillator geometry and the coupling to the phototube must be carefully considered: If the light can travel a variety of path lengths before being collected, an additional contribution to the time resolution will result. Light collection is widely discussed in the literature (refs. 15, 17–19, 28, 29).

Table 4.2 lists the characteristics of several types of photomultipliers. It is noteworthy to observe that the gain of these tubes ranges from  $\sim 0.5 \times 10^6$  to  $2.5 \times 10^6$ . This gain is strongly affected by the age of the tube and the temperature and will change by a factor of  $\sim 2$  for each 100-V change in high voltage.

For response calculations there are several approximations that will aid in a quick "ballpark" answer:

1. Conversion for absorbed energy (eV) to photons (p) for anthracene is  $\sim 70$  eV/p (ref. 35).
2. Conversion of photons to photoelectrons for S-11 photocathode material is  $\sim 10\%$  (ref. 37).
3. 100% of photons are collected on the photocathode and 100% of cathode-emitted photoelectrons are collected on the first dynode.

The function for the total charge output now becomes

$$Q_o = \frac{E}{70} \times P_{c_{eff}} \times \frac{P}{P_a} \times G \times q(\text{coulombs}) \quad (1)$$

where

$Q_o$  = output charge in coulombs

E = absorbed energy in detector in eV

P/P<sub>a</sub> = detector efficiency compared to anthracene from Table 4.1

G = photomultiplier gain from Table 4.2 or from manufacturer's data

q = charge per electron  $1.6 \times 10^{-19}$  coulomb

P<sub>c<sub>eff</sub></sub> = efficiency of photocathode or  $\sim 10\%$  for S-11.

Table 4.2. Photomultiplier Characteristics.

Type	Maximum Outside Dimensions (in.)		Cathode		Dynode System				Max Overall Voltage (V)
	Diameter	Length	Area (cm <sup>2</sup> )	Response	No. of Dynodes	Type	Material	Gain	
<b>RCA</b>									
5819	2-1/4	5-13/16	14.2	S-11	10	Circ. focused	Cs <sub>3</sub> Sb	2.3 x 10 <sup>6</sup>	1250
6342	2-1/4	5-13/16	14.2	S-11	10	Circ. focused	AgMg	0.55 x 10 <sup>6</sup>	1500
6655	2-1/4	5-13/16	14.2	S-11	10	Circ. focused	Cs <sub>3</sub> Sb	2.3 x 10 <sup>6</sup>	1250
6199	1-9/16	4-9/16	7.75	S-11	10	Circ. focused	Cs <sub>3</sub> Sb	2.8 x 10 <sup>6</sup>	1250
6372	2-9/16	7-3/4	80	S-11	10	Circ. focused	Cs <sub>3</sub> Sb	2.5 x 10 <sup>6</sup>	1200
8055	5-1/2	6-3/4	97	S-11	10	Venetian blind	CuBe	1.5 x 10 <sup>6</sup>	2000
<b>Du Mont</b>									
6291	1-1/2	4-1/4	6.4	S-11	10	Box	AgMg	2 x 10 <sup>6</sup>	2100
6292	2-1/16	5-5/8	13.4	S-11	10	Box	AgMg	2 x 10 <sup>6</sup>	2100
6363		6-1/8	31.4	S-11	10	Box	AgMg	2 x 10 <sup>6</sup>	1800
6364	5-1/4	7-1/2	88.8	S-11	10	Box	AgMg	2 x 10 <sup>6</sup>	1800
<b>EMI/US</b>									
9536B	2	4-3/4	14.6	S-11	10	Venetian blind	Cs <sub>3</sub> Sb	3 x 10 <sup>6</sup>	1700
<b>CBS</b>									
CL1012	1-1/2	4-3/4	8	S-11	10	Venetian blind	AgMg	2 x 10 <sup>6</sup>	1750

If the output charge is integrated on a capacitor as it is on this unit, the capacitor voltage becomes:

$$V_c = \frac{Q_o}{C_i} \quad (2)$$

where

- $V_c$  = capacitor voltage in volts
- $Q_o$  = output charge in coulombs
- $C_i$  = the integrating capacitor in farads

In this unit,  $V_c$  is normally integrated on a 550-pF capacitor and is amplified with a gain of 5 by the internal amplifier.

The preamplifier output rise time will be  $\sim 80$  ns or  $2.2\tau$ , whichever is greater, where  $\tau$  is the major decay time constant of the scintillator (Table 4.1).

If the charge output is driven into a low impedance as a transmission line which is terminated in its characteristic impedance ( $Z_o$ ), then the output peak current will be approximately  $I_o = Q_o/\tau Z_o$ . The rise time will be limited to about 20–50 ns by the PM tube. The voltage output under these conditions is  $I_o Z_o = V_o$ ; Therefore:

$$V_o \approx \frac{Q_o}{\tau Z_o} \quad (3)$$

Example:

**Photomultiplier:** RCA-8055, S-11 photocathode, gain  $1 \times 10^5$  at 1.5 kV.

**Scintillator:** NaI.

**Energy Deposited:** <sup>60</sup>Co 1.33 MeV peak; from Eq. (1)

$$Q_o = \frac{1.33 \times 10^6 \text{ eV}}{70} \times 10\% \\ \times 2.4 \times 1 \times 10^6 \times 1.6 \times 10^{-19}$$

or

$$Q_o \approx 0.7 \times 10^{-10} \text{ coulomb.}$$

This charge integrated on 500 pF and amplified by a gain of 5 as in this unit will be:

$$\frac{0.7 \times 10^{-10} \text{ coulomb}}{500 \times 10^{-12} \text{ F}} \times 5 = 0.7 \text{ V}$$

The rise time from Table 4.1 is  $2.2 \tau_2$ , or 550 ns.

The anode pulse into a 50- $\Omega$  terminated line will yield a current pulse of amplitude,

$$I_o \approx \frac{Q_o}{\tau} = \frac{0.7 \times 10^{-10} \text{ coulomb}}{250 \times 10^{-9} \text{ s}}$$

or

$$I_o \approx \frac{Q_o}{T} = \frac{0.7 \times 10^{-10} \text{ coulomb}}{250 \times 10^{-9} \text{ S}}$$

or

$$I_o \approx 0.28 \text{ mA}$$

and again the rise time will be ~20–50 ns, and the output voltage peak will be  $I_o R_o$  or  $V_{out} = 0.28 \text{ mA} \times 50 \Omega = 14 \text{ mV}$ .

### 4.2. TIMING APPLICATIONS

The different specific applications for the unit are essentially limitless, but since it is designed for timing as well as spectroscopy, two of the most often used coincidence systems are discussed and block diagrams given. From these two block diagrams other applications may be formulated by extension.

#### 4.2.1. TYPICAL FAST-SLOW COINCIDENCE SYSTEM

The block diagram of Fig. 4.1 outlines a fast-slow coincidence system. The words “fast-slow” mean that there are essentially two channels of

information retrieval operating in parallel. The fast channel sets the ultimate resolving time of the coincidence circuitry, while the slow channel selects the pulse height range to be accepted from each detector and by means of a slower coincidence requirement combines this with the fast coincidence requirement to yield information having the criteria of the fast resolving time of the fast channel and the energy selection of the slow channel simultaneously.

In the fast channel the timing amplifiers indicated may be the ORTEC 474 or equivalent or may be a combination of fast amplifier/fast discriminator such as the ORTEC 260 Time Pickoff Unit with the 402A Time Pickoff Control to extract fast timing information.

As shown in Fig. 4.1, the time spectrum from the time to amplitude converter is analyzed as the information channel. Of course, the output of the time-to-amplitude converter could be fed to the single channel analyzer to form an alternate fast coincidence channel, which could then feed the coincidence circuit and allow pulse height analysis of either detector channel as desired.




Fig. 4.1. Typical Fast-Slow Coincidence System Using Scintillators.

Figure 4.2 outlines a simple conventional crossover pick-off coincidence system, which is probably the easiest and most versatile method of doing coincidence when the ultimate resolving time is not required. This method is very easy to use; however,

it results in a  $2\tau$  coincidence resolving time which is theoretically worse by a factor of  $\sim 12$  than may be achieved by leading-edge timing such as was indicated in Fig. 4.1.




Fig. 4.2. Conventional Crossover Pickoff Coincidence System.

#### 4.2.2. GAMMA-GAMMA COINCIDENCE SYSTEM FOR THE HIGH-PURITY GERMANIUM (HPGe) DETECTOR

Figure 4.3 is a block diagram of an experimental setup that is quite versatile for studying decay schemes and transitional levels by means of coincidence between an HPGe detector and a scintillation detector. With this block diagram the experimenter may study either energy information or time information associated with the coincidence events.

#### 4.3. SCINTILLATION SPECTROSCOPY

Scintillation spectroscopy implies the measurement of energy by the direct conversion of energy to light in a scintillator and the detection thereof by the photomultiplier. The system to perform this function is one of the most simple, requiring only the phosphor, the photomultiplier, a base structure for the PM tube, a preamplifier, and a linear amplifier with some type of measuring device such as a multichannel analyzer (Fig. 4.4). With this system one may study directly the energy released in the phosphor by some incident radiation.




Fig. 4.3. Gamma-Gamma Coincidence System Using HPGe Detectors.




Fig. 4.4. Scintillation Spectrometry System.

## 5. MAINTENANCE

The resistor divider string of this unit is composed only of passive components; so the only maintenance to be expected is replacement of components that have failed because of age. Table 5.1 lists the approximate dynode voltages for

comparative purposes. Almost all failures of the dynode string may be isolated by removing the PM tube and making these measurements. Use a voltmeter of  $>20,000 \Omega/V$  for this measurement.

Since the amplifier is a high, open-loop gain (~6000) operational amplifier with feedback via R23 and R24, failure of almost any component will cause the output to go to a dc level of  $\sim\pm 5$  V to  $\pm 24$  V.

Troubleshooting involves a careful analysis of the dc levels from Table 5.2. Replacement of the parts is not critical except Q1–Q4 should be very high gain ( $\beta$  or hfe) low-noise NPN transistors. R19, R23, and R24 should be of precision metal film for low noise and good stability.

C7 should be a high quality dipped mica or silver mica capacitor.

D1 and D2 have no role in the amplifier operation and are provided for circuit protection only.

$R19 \times C7$  product is the decay time constant. The 90% to 10% decay time in seconds should be  $\sim 2.2 \times R19(\Omega) \times C7(F)$ .

Dynamic testing can be done with a pulser connected to the test input. With the normal  $C_{test}$  of

100 pF and C7 of 500 pF the preamplifier output amplitude should be approximately the same amplitude and the same polarity as the pulser input. The decay time should be  $\sim 50$ – $150$   $\mu$ s, depending on the pulser. If a long square wave is used for testing, the decay time should be  $\sim 110$   $\mu$ s, 90% to 10%. The rise time should be 80–100 ns.

### 5.1. FACTORY REPAIR

This instrument can be returned to the ORTEC factory for service and repair at a nominal cost. Our standard procedure for repair ensures the same quality control and checkout that are used to a new instrument. Always contact Customer Service at ORTEC, (865) 482-4411, before sending in an instrument for repair to obtain shipping instructions and so that the required Return Authorization Number can be assigned to the unit. Write this number on the address label and on the package to ensure prompt attention when it reaches the factory.

Table 5.1

HV = +2000 V	
PM Socket	Voltage
1	+330
2	+500
3	+660
4	+830
5	+1000
6	+1170
7	+1350
8	+1520
9	+1690
10	+1840
11	+2000
12	NC
13	0–+330 (Focus control).
14	Ground

Table 5.2. Preamplifier Voltages.

Power Cable	Nominal	Volts	
	Volts	(Min)	(Max)
Pin 7	+24	+23	+25
Pin 6	-24	-23	-25
Preamplifier			
Q1b	0	-0.05	+0.05
Q1e	-0.6	-0.4	-0.8
Q1c	+16.4	+15.0	+17.5
Q2e	-1.2	-0.7	-1.5
Q3b	-0.6	-0.4	-0.8
Q4b	0	-0.1	+0.1
Q4c	+4.0	+2.0	+5.0
Q5e	+17.0	+16.0	+18.0
Q6b	+0.6	0	+1.1
Q6e	0	+0.5	-0.5
Q6c	+20.0	+18.0	+22.0
Q7b	-12.0	-11.0	-13.0
Q7e	-12.6	-11.6	-13.6
Q8e	-12.6	-11.6	-13.6

## 6. BIBLIOGRAPHY

1. P.R. Orman, *Nucl. Instr. Methods* 21(1), 121 (1963).
2. D.L. Wieber and H.W. Lefevre, *IEEE Trans. Nucl. Sci.* NS-13(1), 406 (1966).
3. D.A. Gedcke and W. J. McDonald, *Nucl. Instr. Methods* 55(2), 377 (1967).
4. D. A. Gedcke and W. J. McDonald, *Nucl. Instr. Methods* 58(2), 253 (1968).
5. W.J. McDonald and D.A. Gedcke, "Electronics for Fast Neutron Work," presented at the International Symposium on Nuclear Electronics, Versailles, September 1968.
6. R. Nutt, *IEEE Trans. Nucl. Sci.* NS-14(1), 110 (1967).
7. P.D. Compton, Jr., and W.A. Johnson, *IEEE Trans. Nucl. Sci.* NS-14(1) 116 (1967).
8. M. Bertolaccini et al., *Nucl. Instr. Methods* 51(2), 325 (1967).
9. E. Gatti and V. Svelto, *Nucl. Instr. Methods* 43(1), 248 (1966).
10. S. Donati, E. Gatti, and V. Svelto, *Nucl. Instr. Methods* 46(1), 165 (1967).
11. L.G. Hyman, R.M. Schwarcz, and R.A. Schluter, *Rev. Sci. Instr.* 35(3), 393 (1964).
12. L.G. Hyman, *Rev. Sci. Instr.* 36(2), 193 (1965).
13. C.R. Kerns, *IEEE Trans. Nucl. Sci.* NS-14(1), 449 (1967).
14. M. Cocchi and A Rota, *Nucl. Instr. Methods* 55(2), 365 (1967).
15. G. Bertolini et al., *IEEE Trans. Nucl. Sci.* NS-13(3), 119 (1966).
16. J.A. Miede, E. Ostertag, and A. Coche, *IEEE Trans. Nucl. Sci.* NS-13(3), 127 (1966).
17. A. Schwarzchild, *Nucl. Instr. Methods* 21(1), 1 (1963).
18. G. Present et al., *Nucl. Instr. Methods* 31(1), 71 (1964).
19. W. J. McDonald and D.A. Gedcke, *Nucl. Instr. Methods* 55(1), 1 (1967).
20. R.E. Bell, *Nucl. Instr. Methods* 43(2), 211 (1966).
21. R.L. McGuire and R.C. Palmer, *IEEE Trans. Nucl. Sci.* NS-14(1), 217 (1967).
22. F.J. Lynch *IEEE Trans Nucl. Sci.* NS-15(3), 102 (1968).
23. F.T. Kuchnir and F.J. Lynch, *IEEE Trans. Nucl. Sci.* NS-15(3), 107 (1968).
24. W.R. Wall and K. I. Roulston, *IEEE Trans. Nucl. Sci.* NS-15(3), 153 (1968).
25. J. Kirkbride, E. C. Yates, and D.G. Crandall, *Nucl. Instr. Methods* 52(2), 293 (1967).

26. F. J. Lynch, *IEEE Trans. Nucl. Sci.* NS-13(3), 140 (1966).
27. A. Houdayer, S.K. Mark, and R.E. Bell, *Nucl. Instr. Methods* 59(2), 319 (1968).
28. E. Rosenstingl et al., *Nucl. Instr. Methods* 58(1), 61 (1968).
29. P.K.F. Greider, *Nucl. Instr. Methods* 55(2), 295 (1967).
30. J. Braunsfurth and H.J. Körner, *Nucl. Instr. Methods* 34(2), 202 (1965).
31. M.L. Rousch, M.A. Wilson, and W.F. Hornyak, *Nucl. Instr. Methods* 31(1), 112 (1964), and references contained therein.
32. D. Landis and F.S. Goudling, *Nat. Acad. Sci. – Nat. Res. Council Publ.* 1184, 143 (1963).
33. W. J. McDonald and D. A. Gedcke, "Detection System for a Fast Neutron Time-of-Flight Spectrometer," Nuclear Research Center Report, Physics Department, University of Alberta, Edmonton, Alberta, Canada (1967) (copies available from the authors).
34. W.J. Price, *Nuclear Radiation Detection*, 2d ed, McGraw-Hill, New York, 1964.
35. Yu K. Akimov, "Scintillation Counters in High Energy Physics," Academic Press, New York, 1965.
36. D. Gedcke and C.W. Williams, "High Resolution Time Spectroscopy. 1. Scintillation Detectors," ORTEC Publication August 1968.
37. RCA Photocathode Spectral Response Characteristics Chart P IT-701B.