

**Models 9301 and 9301S
Fast Preamplifiers
Operating and Service Manual**

Advanced Measurement Technology, Inc.

a/k/a/ ORTEC[®], a subsidiary of AMETEK[®], Inc.

WARRANTY

ORTEC* warrants that the items will be delivered free from defects in material or workmanship. ORTEC makes no other warranties, express or implied, and specifically NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

ORTEC's exclusive liability is limited to repairing or replacing at ORTEC's option, items found by ORTEC to be defective in workmanship or materials within one year from the date of delivery. ORTEC's liability on any claim of any kind, including negligence, loss, or damages arising out of, connected with, or from the performance or breach thereof, or from the manufacture, sale, delivery, resale, repair, or use of any item or services covered by this agreement or purchase order, shall in no case exceed the price allocable to the item or service furnished or any part thereof that gives rise to the claim. In the event ORTEC fails to manufacture or deliver items called for in this agreement or purchase order, ORTEC's exclusive liability and buyer's exclusive remedy shall be release of the buyer from the obligation to pay the purchase price. In no event shall ORTEC be liable for special or consequential damages.

Quality Control

Before being approved for shipment, each ORTEC instrument must pass a stringent set of quality control tests designed to expose any flaws in materials or workmanship. Permanent records of these tests are maintained for use in warranty repair and as a source of statistical information for design improvements.

Repair Service

If it becomes necessary to return this instrument for repair, it is essential that Customer Services be contacted in advance of its return so that a Return Authorization Number can be assigned to the unit. Also, ORTEC must be informed, either in writing, by telephone [(865) 482-4411] or by facsimile transmission [(865) 483-2133], of the nature of the fault of the instrument being returned and of the model, serial, and revision ("Rev" on rear panel) numbers. Failure to do so may cause unnecessary delays in getting the unit repaired. The ORTEC standard procedure requires that instruments returned for repair pass the same quality control tests that are used for new-production instruments. Instruments that are returned should be packed so that they will withstand normal transit handling and must be shipped PREPAID via Air Parcel Post or United Parcel Service to the designated ORTEC repair center. The address label and the package should include the Return Authorization Number assigned. Instruments being returned that are damaged in transit due to inadequate packing will be repaired at the sender's expense, and it will be the sender's responsibility to make claim with the shipper. Instruments not in warranty should follow the same procedure and ORTEC will provide a quotation.

Damage in Transit

Shipments should be examined immediately upon receipt for evidence of external or concealed damage. The carrier making delivery should be notified immediately of any such damage, since the carrier is normally liable for damage in shipment. Packing materials, waybills, and other such documentation should be preserved in order to establish claims. After such notification to the carrier, please notify ORTEC of the circumstances so that assistance can be provided in making damage claims and in providing replacement equipment, if necessary.

Copyright © 2002, Advanced Measurement Technology, Inc. All rights reserved.

*ORTEC[®] is a registered trademark of Advanced Measurement Technology, Inc. All other trademarks used herein are the property of their respective owners.

CONTENTS

WARRANTY	ii
SAFETY INSTRUCTIONS AND SYMBOLS	iv
SAFETY WARNINGS AND CLEANING INSTRUCTIONS	v
1. DESCRIPTION	1
2. SPECIFICATIONS	1
2.1. PERFORMANCE	1
2.2. CONNECTORS	1
2.3. ELECTRICAL AND MECHANICAL	1
3. INSTALLATION	2
3.1. CONNECTION TO DETECTOR	2
3.2. CONNECTION OF OUTPUT	2
3.3. INPUT POWER	2
3.4. PHYSICAL LOCATION	2
4. OPERATING INSTRUCTIONS	2
5. MAINTENANCE AND REPAIR	3

SAFETY INSTRUCTIONS AND SYMBOLS

This manual contains up to three levels of safety instructions that must be observed in order to avoid personal injury and/or damage to equipment or other property. These are:

DANGER Indicates a hazard that could result in death or serious bodily harm if the safety instruction is not observed.

WARNING Indicates a hazard that could result in bodily harm if the safety instruction is not observed.

CAUTION Indicates a hazard that could result in property damage if the safety instruction is not observed.

Please read all safety instructions carefully and make sure you understand them fully before attempting to use this product.

In addition, the following symbol may appear on the product:



ATTENTION – Refer to Manual



DANGER – High Voltage

Please read all safety instructions carefully and make sure you understand them fully before attempting to use this product.

SAFETY WARNINGS AND CLEANING INSTRUCTIONS

DANGER Opening the cover of this instrument is likely to expose dangerous voltages. Disconnect the instrument from all voltage sources while it is being opened.

WARNING Using this instrument in a manner not specified by the manufacturer may impair the protection provided by the instrument.


Cleaning Instructions

To clean the instrument exterior:

- Unplug the instrument from the ac power supply.
- Remove loose dust on the outside of the instrument with a lint-free cloth.
- Remove remaining dirt with a lint-free cloth dampened in a general-purpose detergent and water solution. Do not use abrasive cleaners.

CAUTION To prevent moisture inside of the instrument during external cleaning, use only enough liquid to dampen the cloth or applicator.

- Allow the instrument to dry completely before reconnecting it to the power source.



ORTEC 9301 and 9301S FAST PREAMPLIFIERS

1. DESCRIPTION

The ORTEC 9301 and 9301S Fast Preamplifiers are designed for use in fast timing system applications. They feature very low noise and fast risetime to facilitate precise timing measurements. Their small size and light weight permit them to be used in close proximity to the detectors they serve, and this minimizes noise in the output signals. The two Preamplifiers are identical except for the mating power connector at the remote end of the power cable. Although the reference in this manual is to the 9301, the instructions also relate to the 9301S.

The 9301 can be used with photomultipliers, electron multipliers, and other detectors that are employed in photon counting or ion counting techniques. When the preamplifier is connected to the detector, the 50Ω input impedance of the 9301 becomes the load resistance for the detector output current pulse.

The risetime of the 9301 output pulse is <1.5 nsec and the voltage gain is 10. This combination ensures the user that each input pulse from the detector will rise quickly through a normal discriminator threshold for a precise timing identification. The voltage gain of 10 also enhances the signal-to-noise ratio for the system.

Power for operation of the 9301 must be obtained from either the instrument to which its output pulses are furnished or from an ORTEC 114 Preamplifier Power Supply. A 10-ft power cable is built into the unit and terminated with an Amphenol 17-20090 plug. The Amphenol plug is compatible with the Preamp power source on a 114 Preamplifier Power Supply, a 9302 Amplifier/Discriminator, or any ORTEC NIM-standard amplifier module. For the 9301S the power cable is terminated with a DIN plug, compatible with the ORTEC Brookdeal 5C1 Photon Counting System.

2. SPECIFICATIONS

2.1. PERFORMANCE

Input Impedance 50Ω

Voltage Gain 10 ($\pm 2\%$) noninverting.

Risetime <1.5 nsec.

Input RMS Noise Equivalent <25 μ V.

Output Impedance 50Ω nominal.

Output Dynamic Range $>\pm 0.7$ V into 50Ω

Nonlinearity $<1\%$.

Temperature Gain Instability $<0.1\%/^{\circ}\text{C}$.

2.2. CONNECTORS

Input Connector BNC, front panel.

Output Connector BNC, rear panel.

2.3. ELECTRICAL AND MECHANICAL

Power Required

+12 V, 30 mA;

-12 V, 30 mA.

Specify the 9301S for use with 5C1 Photon Counting System.

Specify the 9301 for connection to all other ORTEC equipment.

Dimensions 1.25 X 2 X 2.875 in. plus 10-ft cable.

3. INSTALLATION

3.1. CONNECTION TO DETECTOR

Connect an ac-coupled output from the detector to the Input BNC on the 9301 with coaxial cable. This cable should be kept as short as practical to minimize the physical distance through which the detector output current pulse has to pass. The amount of stray noise that will be added to the signal increases as the cable length increases. Use standard 50 Ω coaxial cable.

An input pulse to the 9301 can be of either polarity; so it can originate from the anode or any dynode in, for example, a photomultiplier tube circuit. It must be ac-coupled to prevent high voltage from being applied across the 50 Ω dc-coupled input circuit of the preamplifier. The output pulse from the 9301 will have the same polarity as the input and will be amplified by a factor of 10. The linear range of the output is +700 mV to -700 mV.

3.2. CONNECTION OF OUTPUT

The output impedance of the 9301 is nominally 50 Ω . It is designed to drive any counter, amplifier, or timing instrument through a long length of 50 Ω cable. The input impedance of the driven instrument is usually on the order of 1000 Ω or more, and the series termination provided by the 9301 will usually be sufficient to prevent pulse reflections. If interference results from the basic interconnection, it can be reduced by an additional shunt termination of the 50 Ω at the input to the driven instrument, but this will also attenuate the signal amplitude by 50%.

Since the power cable for the 9301 is 10 ft. long (standard), the separation of the preamplifier and the instrument being driven will usually be 10 ft maximum. Do not use signal cables that are excessive in length; the signal cable will normally be no longer than 10 ft.

3.3. INPUT POWER

Power for the 9301 is supplied through the captive 10-ft power cable. It can be furnished from any of a variety of sources, such as the ORTEC instrument to which the output pulses are being furnished or an ORTEC 114 Preamplifier Power Supply. The power requirements are +12 V, 30 mA, and -12 V, 30 mA, and can be furnished from a center-tapped 24-V battery if desired.

The 9301 Fast Preamplifier includes an Amphenol 17-20090 plug, at the remote end of the power cable, which is wired to accept +12 V through pin 4, -12 V through pin 9, and ground through pins 1 and 2. It is directly compatible with any ORTEC NIM-standard module that has a Preamp power connector mounted on the rear panel, or an ORTEC 114 Preamplifier Power Supply, a 9302 Amplifier/Discriminator, or a 9305 Photometer.

The 9301S Fast Preamplifier includes a DIN plug at the remote end of the power cable. This plug is compatible with the preamplifier power source that is included in the ORTEC Brookdeal 5C1 Counting System.

Power cable lengths other than the standard 10 ft are available on special order, with either an Amphenol 17-20090 or a DIN plug used for termination.

3.4. PHYSICAL LOCATION

The 9301 weighs only 6 oz and can normally be supported by its cables. There are no mounting facilities required for the unit, and it will usually be located in close proximity to the source of input pulses to minimize stray noise pickup.

4. OPERATING INSTRUCTIONS

There are no operating controls or adjustments for the 9301 Fast Preamplifier. Simply connect its power and input and output signals as described in Section 3 and its operation will be automatic.

If desired, a pulse generator can be substituted for the source of input pulses to check system response. The input pulse amplitude to the 9301 should not exceed 70 mV of either polarity.

5. MAINTENANCE AND REPAIR

Keep the signal and power connectors clean and free from any accumulation of dust. No other maintenance should be required. If the preamplifier is suspected of malfunction, check to see that the

required dc power levels are being furnished and then use a pulse generator to check its response to input signals as suggested in Section 4.

