

**Model VT120
Fast Timing Preamplifier
Operating and Service Manual**

Advanced Measurement Technology, Inc.

a/k/a/ ORTEC[®], a subsidiary of AMETEK[®], Inc.

WARRANTY

ORTEC* warrants that the items will be delivered free from defects in material or workmanship. ORTEC makes no other warranties, express or implied, and specifically NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

ORTEC's exclusive liability is limited to repairing or replacing at ORTEC's option, items found by ORTEC to be defective in workmanship or materials within one year from the date of delivery. ORTEC's liability on any claim of any kind, including negligence, loss, or damages arising out of, connected with, or from the performance or breach thereof, or from the manufacture, sale, delivery, resale, repair, or use of any item or services covered by this agreement or purchase order, shall in no case exceed the price allocable to the item or service furnished or any part thereof that gives rise to the claim. In the event ORTEC fails to manufacture or deliver items called for in this agreement or purchase order, ORTEC's exclusive liability and buyer's exclusive remedy shall be release of the buyer from the obligation to pay the purchase price. In no event shall ORTEC be liable for special or consequential damages.

Quality Control

Before being approved for shipment, each ORTEC instrument must pass a stringent set of quality control tests designed to expose any flaws in materials or workmanship. Permanent records of these tests are maintained for use in warranty repair and as a source of statistical information for design improvements.

Repair Service

If it becomes necessary to return this instrument for repair, it is essential that Customer Services be contacted in advance of its return so that a Return Authorization Number can be assigned to the unit. Also, ORTEC must be informed, either in writing, by telephone [(865) 482-4411] or by facsimile transmission [(865) 483-2133], of the nature of the fault of the instrument being returned and of the model, serial, and revision ("Rev" on rear panel) numbers. Failure to do so may cause unnecessary delays in getting the unit repaired. The ORTEC standard procedure requires that instruments returned for repair pass the same quality control tests that are used for new-production instruments. Instruments that are returned should be packed so that they will withstand normal transit handling and must be shipped PREPAID via Air Parcel Post or United Parcel Service to the designated ORTEC repair center. The address label and the package should include the Return Authorization Number assigned. Instruments being returned that are damaged in transit due to inadequate packing will be repaired at the sender's expense, and it will be the sender's responsibility to make claim with the shipper. Instruments not in warranty should follow the same procedure and ORTEC will provide a quotation.

Damage in Transit

Shipments should be examined immediately upon receipt for evidence of external or concealed damage. The carrier making delivery should be notified immediately of any such damage, since the carrier is normally liable for damage in shipment. Packing materials, waybills, and other such documentation should be preserved in order to establish claims. After such notification to the carrier, please notify ORTEC of the circumstances so that assistance can be provided in making damage claims and in providing replacement equipment, if necessary.

Copyright © 2002, Advanced Measurement Technology, Inc. All rights reserved.

*ORTEC[®] is a registered trademark of Advanced Measurement Technology, Inc. All other trademarks used herein are the property of their respective owners.

CONTENTS

| | |
|-------------------------------------------------|----|
| WARRANTY | ii |
| SAFETY INSTRUCTIONS AND SYMBOLS | iv |
| SAFETY WARNINGS AND CLEANING INSTRUCTIONS | v |
| 1. DESCRIPTION | 1 |
| 2. SPECIFICATIONS | 1 |
| 2.1. PERFORMANCE | 1 |
| 2.2. ELECTRICAL AND MECHANICAL | 1 |
| 2.3. ACCESSORIES | 1 |
| 3. INSTALLATION | 2 |
| 3.1. GENERAL | 2 |
| 4. CIRCUIT DESCRIPTION | 2 |
| 5. MAINTENANCE | 2 |

SAFETY INSTRUCTIONS AND SYMBOLS

This manual contains up to three levels of safety instructions that must be observed in order to avoid personal injury and/or damage to equipment or other property. These are:

- DANGER** Indicates a hazard that could result in death or serious bodily harm if the safety instruction is not observed.
- WARNING** Indicates a hazard that could result in bodily harm if the safety instruction is not observed.
- CAUTION** Indicates a hazard that could result in property damage if the safety instruction is not observed.

Please read all safety instructions carefully and make sure you understand them fully before attempting to use this product.

In addition, the following symbol may appear on the product:



ATTENTION—Refer to Manual



DANGER—High Voltage

Please read all safety instructions carefully and make sure you understand them fully before attempting to use this product.

SAFETY WARNINGS AND CLEANING INSTRUCTIONS

DANGER Opening the cover of this instrument is likely to expose dangerous voltages. Disconnect the instrument from all voltage sources while it is being opened.

WARNING Using this instrument in a manner not specified by the manufacturer may impair the protection provided by the instrument.


Cleaning Instructions

To clean the instrument exterior:

- Unplug the instrument from the ac power supply.
- Remove loose dust on the outside of the instrument with a lint-free cloth.
- Remove remaining dirt with a lint-free cloth dampened in a general-purpose detergent and water solution. Do not use abrasive cleaners.

CAUTION To prevent moisture inside of the instrument during external cleaning, use only enough liquid to dampen the cloth or applicator.

- Allow the instrument to dry completely before reconnecting it to the power source.



ORTEC MODEL VT120 FAST TIMING PREAMPLIFIER

1. DESCRIPTION

The VT120 preamplifier is a high-performance, wide-bandwidth preamplifier designed for boosting very fast linear signals from photomultipliers, electron multipliers, silicon surface-barrier detectors, and other detectors used in fast timing applications. The rise time on all versions is <1 ns with a 5-V output, enabling excellent timing resolution.

The VT120 is a single-channel unit in a small preamp package. It is available with a gain of 200, noninverting (A version); a gain of 200, inverting (B version); or a gain of 20, noninverting (C version). BNC connectors are used for signal connections on the VT120. A cable is available (Model C-VT120) for connecting power between the Model VT120 and conventional preamplifier power outputs using Amphenol-type 17-80090 connectors.

2. SPECIFICATIONS

2.1. PERFORMANCE

GAIN (10% gain tolerance on all versions):

A Version 200, noninverting.

B Version 200, inverting.

C Version 20, noninverting.

RISE TIME ≤ 1 ns.

NOISE ≤ 20 μ V rms equivalent input noise.

BANDWIDTH 10 to 350 MHz.

OUTPUT RANGE 0 to -5 V with 50- Ω load.

INPUT BNC connector; input impedance 50 Ω .

OUTPUT BNC connector; 0 to -5 V output with a 50- Ω load. Output impedance ≤ 1 n.

2.2. ELECTRICAL AND MECHANICAL

POWER REQUIRED +12 V, 50 mA (uses LEMO power connector that is compatible with accessory cable C-VT120).

DIMENSIONS Aluminum housing 5.8 X 5.1 X 1.6 cm (2.3 X 2.0 X 0.63 in.).

WEIGHT 0.2 kg (0.4 lb).

2.3. ACCESSORIES

C-VT120 cable assembly with connections between VT120 power input (LEMO) and Amphenol-type 17-80090 preamplifier power connectors that are compatible with other ORTEC NIM-standard modules.

3. INSTALLATION

3.1. GENERAL

The VT120 contains no internal power supply and must obtain power from a NIM-standard bin and power supply, usually by connecting the C-VT120 accessory power cable between the VT120 and a NIM-standard module with a mating Amphenol-type 17-80090-15 preamp power connector. The bin and

power supply should be turned Off when modules are inserted or removed. The power supply voltages should be checked after modules are inserted. Ensure that the VT120 has sufficient cooling air circulating to prevent any localized heating of the solid-state circuitry used throughout the unit. The VT120 should not be subjected to temperatures in excess of 50°C.

4. CIRCUIT DESCRIPTION

The VT120 preamplifier is available in three versions. The "A" version has a noninverting gain of 200, the "B" version has an inverting gain of 200, and the "C" version has a noninverting gain of 20. The schematic diagram for each version is attached at the end of this manual.

The schematic diagram for the VT120A at the end of this manual shows that the "A" version has five gain stages, Q1 through Q4, connected in the Common-Emitter configuration and Q5 connected in the Emitter-Follower configuration. Each stage is ac-coupled to isolate bias voltages and to strongly reject low-frequency components. Input protection is furnished by diode D1.

Consider the first stage formed around Q1. Resistor R2 furnishes base-injection bias for the transistor, the series combination C2 and R3 provides shunt feedback, and resistor R4 is an unbypassed emitter resistor which, along with C3, peaks the high-frequency response of the stage. Each of the other Common-Emitter stages operates in a similar manner with specific component values set for optimum gain distribution and rise time. The final two Common-Emitter stages, Q3 and Q4, have variable peaking capacitors, C20 and C21, which

are set in test to give an overall rise time of <1 ns. Transistor Q5 is connected in the Emitter-Follower configuration and provides excellent output drive capability, 0 to -5 V, and low output impedance, <1 Ω.


The schematic diagram of the VT120B is attached at the end of this manual. The "B" version has one less Common-Emitter amplifier stage than does the "A" version. This causes an overall signal inversion from input to output. The operation of each stage is similar to that of the "A", version except for the distribution of gain and rise time among the various preamplifier stages. Peaking capacitors C16 and C17 are set in test to give an overall rise time <1 ns.

The schematic diagram of the VT120C is also attached at the end of this manual. The "C" version has two less Common-Emitter amplifier stages than does the "A" version. There is no signal inversion from input to output, and the overall gain is nominally 20. The operation of each stage is similar to that of the "A" version except for the distribution of gain and rise time among the various preamplifier stages. Peaking capacitors C3 and C7 are set in test to give an overall rise time <1 ns.

5. MAINTENANCE


The VT120 Fast Timing Preamplifier requires very little maintenance other than routine removal of dust and tightening of mechanical connections.

1. RESISTOR VALUES ARE IN OHMS.
2. RESISTOR R59 IS SELECTED DURING TEST.



| | | |
|-------|-----|----------------------------------------------|
| DATE | REV | DESCRIPTION |
| 12/27 | 1 | ORTEC |
| 12/27 | 2 | 100 MIDLAND ROAD, OAK RIDGE, TENNESSEE 37830 |
| | | MODEL VT1250A |
| | | 031A0, SCHENK, FAST FINDING AMP. |
| | | DESIGNED BY: J. J. HALL |
| | | DESIGNED DATE: 12/27/50 |
| | | TEST NUMBER: 747093 |
| | | TEST DATE: 12/27/50 |
| | | TESTER: J. J. HALL |
| | | TEST ROOM: 747093 |
| | | SHEET 1 OF 1 |

- 1. RESISTOR VALUES ARE IN OHMS.
- 2. R6 IS SELECTED DURING TEST.



| | | | |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|-----|--------|
| DATE | REV | BY | DR/TEC |
| 7/77 | 01/21 | 100 | 100 |
| MODEL VT120C 0140, SCHEM. FAST TIMING AMP. FOR USE WITH THE PERFORMANCE OF THE 100 SERIES OF THE 75000 SERIES OF THE 75000 T. J. ALLAN 747100 747100 | | | |